



始 转录本拼装

为什么转录本可以拼接起来?








有参组装和无参组装








有参转录组分析流程 - 可变剪接








选择性剪接丰富了蛋白产物的多样性







卵子特有的Dicer剪切异构体控制RNAi通路



ARTICLE | VOLUME 155, ISSUE 4, P807-816, NOVEMBER 07, 2013

A Retrotransposon-Driven Dicer Isoform Directs Endogenous Small Interfering RNA Production in Mouse Oocytes

Matyas Flemr • Radek Malik • Vedran Franke • ... Radislav Sedlacek • Kristian Vlahovicek •

Highlights

- Retrotransposition gave rise to an oocyte-specific Dicer isoform (Dicer^O) in mice
- Dicer^O is N-terminally truncated and has higher activity than somatic Dicer
- Dicer^O controls the endogenous RNAi pathway and is essential for mouse oocytes
- Low Dicer activity and low dsRNA abundance constrain endogenous RNAi in mammals
 易生信, 毕生缘: 培训版权所有。




二代测序难获得准确的转录本




- 二代测序难拼出准确转录本,更适合检验单个剪接 位点的变化






Transcript classification codes









可变剪接分析

可变剪接的类型




Alternative Splicing Events

Skipped exon (SE)




外显子跳跃

Alternative 5' splice site (A5SS)



可变的5' 剪接位点

Alternative 3' splice site (A3SS)



可变的3' 剪接位点

Mutually exclusive exons (MXE)




外显子互斥

Retained intron (RI)



内含子保留








Alternatively spliced exon

样品特异的选择性剪接事件的鉴定







rMATS必需参数



Required Parameters:

--bi STARIndexFolder

--s1 s1.txt A text file contains FASTQ file(s) for the sample_1.(Only if using fastq)

--s2 s2.txt A text file contains FASTQ file(s) for the sample_2.(**Only if using** fastq)

--b1 b1.txt A text file records mapping results for the sample_1 in bam format.

(Only if using bam)

--b2 b2.txt A text file records mapping results for the sample_2 in bam format. (Only if using bam)

-t readType

Type of read used in the analysis. readType is either 'paired' or 'single'.

'paired' is for paired-end data and 'single' is for single-end data

--readLength <int> The length of each read

--gtf gtfFile An annotation of genes and transcripts in GTF format

The folder name of the STAR binary indexes (i.e., the name of the folder that contains SA file). For example, use ~/STARindex/hg19 for hg19. Only if using fastq)

--od outDir The output directory

s1.txt和s2.txt是以逗号分隔的样本FASTQ文件

b1.txt和b2.txt是以逗号分隔的样本bam文件

-t readType 双端测序为 paired , 单端测序为 single

测序reads的长度

gtf注释文件

输入文件是fastq格式时, 指定STAR索引文件位置

输出结果文件夹



rMATS可选参数



\sim	-	41	_	-	-	١.
О	D	u	o	п	a	ı.
_	•	-	_		_	•

	The "anchor length" or "overhang length" used in the aligner. At least
tophatAnchor <int></int>	"anchor length" NT must be mapped to each end of a given junction.
	The default is 6. (This parameter applies only if using fastq)

- --nthread <float> The number of thread. The optimal number of thread should be equal to the number of CPU core.
- The cutoff splicing difference. The cutoff used in the null hypothesis test for differential splicing. The default is 0.0001 for 0.01% difference. Valid: 0 ≤ cutoff < 1
- --tstat <float> The number of thread for statistical model.
- --statoff Turn statistics part off.
- Library type. Default is unstranded (fr-unstranded). Use fr-firststrand

建议加上

线程数

统计计算

时线程数



-libType libraryType

or fr-secondstrand for strand-specific data.

rMATS运行脚本



- o python rMATS-turbo-xxx-UCSx/rmats.py --s1 s1.txt --s2 s2.txt --gtf gtfFile --bi STARindexFile --od out_directory -t paired --nthread 2 --tstat 2 --readLength 101 --tophatAnchor 8 --cstat 0.0001 #输入文件是fastq格式
- o python rMATS-turbo-xxx-UCSx/rmats.py --b1 b1.txt --b2 b2.txt -gtf gtfFile --od bam_out_directory -t paired --nthread 2 --tstat 2 -- readLength 101 --cstat 0.0001 --libType fr-unstranded #輸入文件是bam格式



rMATS输出结果解释



Table 11.1: rMATS输出结果解释

V1	V2		
AS_Event.MATS.JC.txt	依赖Junction reads定义的剪接位点丰度		
AS_Event.MATS.JCEC.txt	依赖Junction reads和对应原件其它reads定义的剪接位点丰度		
fromGTF.AS_Event.txt	从GTF中提取的所有的可能的AS(alternative splicing)事件		
JC.raw.input.AS_Event.txt	剪接位点的Junction reads count		
JCEC.raw.input.AS_Event.txt	剪接位点和剪接元件上的总reads count		



rMATS运行结果




- AS_Event.MATS.JC.txt
 evaluates splicing with only reads that span splicing junctions
- AS_Event.MATS.JCEC.txt
 evaluates splicing with reads that span splicing junctions and reads on target
- fromGTF.AS_Event.txt
 all possible alternative splicing (AS) events derived from GTF and RNA
- JC.raw.input.AS_Event.txt
 evaluates splicing with only reads that span splicing junctions
- JCEC.raw.input.AS_Event.txt
 evaluates splicing with reads that span splicing junctions and reads on target



rMATS输出结果每列含义解







rMATS输出结果每列含义解释



ID	可变剪接分析编号,标示使用,无特殊含义		1stExonStart 0base	MXE中第一个使用的外显子的起始位置
GeneID	注释的基因的名字		1stExonEnd	MXE中第一个使用的外显子的终止位置
geneSymbol	基因名字			
chr	染色体名字		2ndExonStart_0base	MXE中第二个使用的外显子的起始位置
strand	基因处在染色体的正链还是负链		2ndExonEnd	MXE中第二个使用的外显子的终止位置
riExonStart 0base	保留内含子后的新外显子在基因组的起始位置 (0-based, included)		IJC_T40	T40样品中包含可变剪接位置处Junction reads数,不同的重复用逗号分开
_			SJC_T40	T40样品中去除可变剪接位置处Junction reads数,不同的重复用逗号分开
riExonEnd	保留内含子后的新外显子在基因组的终止位置 (0-based, not included)		IJC_T45	T45样品中包含可变剪接位置处Junction reads数,不同的重复用逗号分开
upstreamES	发生可变剪接的外显子上游外显子起始位置 (0-based, included) 发生可变剪接的外显子上游外显子终止位置 (0-based, not included)		SJC_T45	T45样品中去除可变剪接位置处Junction reads数,不同的重复用逗号分开
upstreamEE			IncFormLen	包含可变剪接部分的转录本的长度
downstreamES	发生可变剪接的外显子下游外显子起始位置 (0-based, included)			
downstreamEE	发生可变剪接的外显子下游外显子终止位置 (0-based, not included)		SkipFormLen	不包含可变剪接部分的转录本长度
longExonStart_0base	A3SS中下游外显子最左侧剪接位点或A5SS中上游外显子起始		Pvalue	差异检验显著性度量
longExonEnd	A3SS中下游外显子终止位置或A5SS中上游外显子最右侧剪接位点		FDR	多重假设检验校正的假阳性率
	A3SS中下游外显子次左侧剪接位点或A5SS中上游外显子起始 A3SS中下游外显子终止位置或A5SS中上游外显子次右侧剪接位点		IncLevel1	样品1 (T40)中包含可变剪接部分的转录本的表达,不同的重复用逗号分开
shortES			IncLevel2	样品2 (T45)中包含可变剪接部分的转录本的表达,不同的重复用逗号分开
shortEE			IncLevelDifference	样品1(T40)的转录本表达量平均值减去样品2的转录本表达平均值(T45)
flankingES	A3SS中上游外显子起始位置或A5SS中下游外显子起始位置			, , , , , , , , , , , , , , , , , , , ,
flankingEE	A3SS中上游外显子终止位置或A5SS中下游外显子终止位置		Remaining columns 占 いいか以作人という 行。	Function annotation

rmats2sashimiplot必需参数



Required Parameters

--s1 s1_rep1.sam[,s1_rep2.sam]

Mapping results for the sample_1 in sam format. Replicates must be in a comma separated list (Only if using sam).

--s2 s2.rep1.sam[,s2.rep2.sam]

Mapping results for the sample_2 in sam format. Replicates must be in a comma separated list (Only if using sam).

--b1 s1_rep1.bam[,s1_rep2.bam]

Mapping results for the sample_1 in bam format.
Replicates must be in a comma separated list
(Only if using bam).

--b2 s2.rep1.bam[,s2.rep2.bam]

Mapping results for the sample_2 in bam format. Replicates must be in a comma separated list (Only if using bam).

-t eventType

Type of event from rMATS result used in the eventType is 'SE', 'A5SS', 'A3SS', 'MXE' or 'RI'. 'SE' is for skipped exon events, 'A5SS' is for alternative 5' splice site events, 'A3SS' is for alternative 3' splice site events, 'MXE' is for mutually exclusive exons events and 'RI' is for retained intron events (Only if using rMATS format result as event file).

-e eventsFile

The rMATS output event file (Only if using rMATS format result as event file).

-c coordinate:annotaionFile

The coordinate of genome region and an annotation of genes and transcripts in GFF3 format. Coordinate and annotation file must be colon separated (Only if using coordinate and annotation file). The label for first sample.

--11 SampleLabel1

--12 SampleLabel2 The label for second sample.

-o outDir

The output directory.

--s1和--s2参数: 逗号分隔的样本 sam文件

--b1和--b2参数: 逗号分隔的样本 bam文件

-t 可变剪接类型

-e 对应剪接类型的 rMATS输出结果

--l1 和--l2:第一组和 第二组样本名称



Rmats2sashimiplot可选参数



Optional:

--exon s <int> The size of scale down exons. The default is 1.

--intron s <int> The size of scale down introns. For example, if

-intron s is 5, it means the size of intron is 5:1

(if the real size of intron is 5, the size in the

plot will be scaled down to 1). The default is 1.

--group-info If user want to divide samples into groups,

they can specify this parameter with a "*.gf" file.

Format specification can be found in following

section.

--min-counts If the junction count is smaller(<) than this float

> number, then this junction would be omitted. The default value is 3. If you want to display all the

numbers, then set it as 0.

--color User can customerize the colors of the plot using a

sequence of color. The number of the colors are

supposed to be corresponding to that of bam files.

eg: --color #FFCC99, #99CC99, #99CC99

--font-size Change the default font size which equals to 8.

--no-text-background Transparent text background.

Hide the numbers of junction.

--hide-number

Print this help message and exit(also --help).

rmats2sashimiplot运行脚本



o rmats2sashimiplot --s1 sampleA.R1.sam,sampleA.R2.sam,sampleA.R3.sam --s2 sampleB.R1.sam,sampleB.R2.sam,sampleB.R3.sam -t SE -e/MATS_output/sampleA_sampleB.SE.MATS.events.txt --l1 sampleA --l2 sampleB --exon_s 1 --intron_s 5 -o events_output #輸入文件是sam格式



rmats2sashimiplot运行脚本




o rmats2sashimiplot --b1 sampleA.R1.bam,sampleA.R2.bam,sampleA.R3.bam --b2 sampleB.R1.bam,sampleB.R2.bam,sampleB.R3.bam -c chr16:-:24944500:24955500:.gff3File --l1 sampleA --l2 sampleB --exon_s 1 --intron_s 5 -o events_output #输入文件是bam格式,用户提供坐标和gff3格式注释文件



sashimiplot











以组的方式展示结果




rmats2sashimiplot --b1 sampleA.R1.bam,sampleA.R2.bam,sampleA.R3.bam
 --b2 sampleB.R1.bam,sampleB.R2.bam,sampleB.R3.bam
 -t SE -e
 /MATS_output/sampleA_sampleB.SE.MATS.events.txt --I1 sampleA --I2 sampleB --exon_s 1 --intron_s 5 -o events_output --group-info grouping.gf

#grouping.gf:

group1name: 1-3

group2name: 4-6



chr16:24953308:24953464:-@chr16:24950685:24950918:-@chr16:24946791:24946960:-



..., 毕生缘; 培训版权所有。

Sequencing costs a lot and gains more







